

## *Certificate of Analysis*

### 1. **Manufacturer**

Mercodia AB, Sylveniusgatan 8A, 754 50 Uppsala, SWEDEN

### 2. **Description**

Catalog no: 10-1195-01

Product: Diabetes Sample Buffer

Lot no: 36341

Expiry date: 2027-10-08

### 3. **Quality control**

Quality control has been performed for lot no 36341 according to standard operating procedures at Mercodia AB, and the product is released based on fulfillment of established acceptance criteria.

### 4. **Intended use**

Mercodia Diabetes Sample Buffer is used for sample dilution for the following Mercodia ELISA assays: Human Insulin, Human Proinsulin, Iso-Insulin, Rat and Mouse Insulin, Rat/Mouse Proinsulin, Rat C-peptide and Animal Insulin.

### 5. **Storage and handling**

Recommended storage of Diabetes Sample Buffer is 2-8°C.

### 6. **Hazardous information**

Please refer to the Material Safety Data Sheet for hazard identification.

### 7. **Quality standard documentation**

The Mercodia Quality System fulfills the QSR and the requirements for the European CE mark, the Directive 98/79/CE on *in vitro* diagnostic medical devices. Mercodia is ISO 13485 certified, for compliance with the International Quality Management System for medical devices.

### 8. **Names and signatures of certifying officers**

Date of analysis:  
2024-10-18.....

Performed by:  
Matias Jäder gm. Jonan Klide 

Date of approval:  
2024-11-28.....

Approved by:  
Elin Westberg EW 