

Certificate

Stability Testing of Mercodia's ELISA Assays


To Whom it May Concern:

Mercodia refers to the below stability testing of their ELISA kits;

Stressed stability studies have been performed on the Mercodia assays where the specifications are still met after up to 8 weeks in +30 °C.

Transport stability has also been studied; After 4 days in +30 °C, 2 days in room temperature and finally +2-8 °C until expiry date, the specifications are still met.

Uppsala November 15, 2021



Pelle Jadeborg
Managing Director Quality
Mercodia AB